

Aesculap® Recall Set

Titanium cures for dental implant cleaning



Aesculap Dental

Recall Set for Dental Implant Cleaning

To ensure the long-term success of a dental implant, consistent aftercare is very important.

To eliminate plaque on the implant-abutment interface without any infection risk, titanium curettes are the optimal choice.



Titanium curettes for dental implant cleaning:

- Very delicate working tips
- Can be re-sharpened
- Biocompatible



DB180T LANGER Titanium curette, fig. 1/2, for premolars, molars, universal; mandible

1/4



DB181T LANGER Titanium curette, fig. 3/4, for premolars, molars, universal; maxilla

1/4



DB182T LANGER Titanium curette, fig. 5/6, for anterior teeth, universal

1/4



Force-controlled periodontal probes


- Defined probing force of 0.2N (20g)
- Titanium Nitride coating
- Exact marking and smoothing



DB764R UNC 15, marking 11111-11111-11111mm



DB765R WHO; PCP 11,5, marking 3,5 – 2 – 3 – 3 mm



DA090 Blue

1/1



DA091 Grey

1/1



DA092 Green

1/1



DA093 Black

1/1

DA090 -DA093 Plastic mouth mirror handles, sterilizable to 240 °C, 135 mm, 5 1/4"



DA074R Stainless steel mouth mirror handle, 140 mm, 5 1/2"

1/1



DA076R Stainless steel mouth mirror handle, 35 mm, 5 1/4"



1/1



DA036R RHODIUM mouth mirror, Ø 22 mm, sales unit: **PAK** = pack of 12 pieces


DA037R RHODIUM mouth mirror, Ø 24 mm, sales unit: **PAK** = pack of 12 pieces

1/1




Aesculap Dental trays

- Safe & effective
- Practical & economical




JG387R Tray, 188 x 84 x 30 mm,
incl. 2 silicone instrument racks (JG394), for 5 instruments

JG388R Tray, 188 x 125 x 30 mm,
incl. 2 silicone instrument racks (JG395), for 9 instruments



JG625R Holding device for identification labels

ICD-002 Plastic identification labels, white,
with 10-digit barcode for praxiDoc or DIOS software applications,
labeling as specified



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Tel. +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | www.aesculap-dental.com

Aesculap – a B. Braun company

'Aesculap' is a registered mark of Aesculap AG.
Subject to technical changes. All rights reserved.
This brochure may only be used for the exclusive
purpose of obtaining information about our
products. Reproduction in any form partial or
otherwise is not permitted.

Brochure no. D12202

0311/3.0/1